



uSpO₂ (USB)

Pulse Oximetry Cable



LNCS® Compatible Connector
(FRONT)

Cable	Typical
Diameter	0.185 inch ± 0.004 inch
Conductors	5
Overall length of cable	6.0 feet max
Color	White

Connector Dimensions	Typical
Width	1.20 inch
Depth	3.60 inch
Thickness	0.90 inch

Power Specifications		
	Typical	Max
Input Voltage	5.0 V	±5%
Average Power	125 mW	250 mW
Line / load regulation	-	50 mV rms

	Connector Description	Number of Pins
Sensor Connector	LNCS Compatible, with Mechanical Latch	9
USB Connector	USB Type A Plug	4

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

